

**Jadavpur University****1<sup>st</sup> BSc 1<sup>st</sup> Sem GE Practical 2019 (new)****Paper: GE/Chem/PR/01****Full Marks: 50**

Q1. Estimate the amount of Fe(II) and Fe(III) present together in a mixture.	30
Q2. Viva-Voce	10
Q3. Laboratory Notebook	10

## Hints Q1

(a) Transfer quantitatively the given sample marked XX into a 250 mL volumetric flask and make up the volume with distilled water to prepare the stock solution.

## (b) (i) Determination of Fe(II)

Take an aliquot of 25 mL of the stock solution in a 250 mL conical flask, dilute with an equal volume of 4(N) sulphuric acid and add 5 mL syrupy phosphoric acid and 2-3 drops of BDS indicator solution. Titrate with supplied standard (N/20)  $K_2Cr_2O_7$  solution, while swirling the flask gently, until the solution changes from bright green to intense red-violet. Record the titre value of dichromate solution.

## (ii) Determination of Fe(II) + Fe(III)

Take an aliquot of 25 mL of the stock solution in a 500 mL conical flask. Add 25 mL conc. HCl. The solution assumes a reddish brown colour. Heat the solution to nearly boiling and then reduce the solution with  $SnCl_2$  solution adding dropwise with constant shaking until the yellow colour of the solution just discharges. Add one drop  $SnCl_2$  in excess. Cool the flask rapidly under tap to room temperature. Add 10 mL 5%  $HgCl_2$  solution at a time, shake and allow to stand for about 5 min when a silky white ppt appears. Dilute the solution with 100 mL of distilled water, add 5 mL syrupy phosphoric acid and 2-3 drops of BDS indicator solution. Titrate with supplied standard (N/20)  $K_2Cr_2O_7$  solution, while swirling the flask gently, until the solution changes from bright green to intense red-violet. Record the titre value of dichromate solution.

(iii) Report the amount of Fe(II) and Fe(III) in the mixture separately.